Product Specification Sheet

Photodiodes



1 motodroux

Description

The Thorlabs FGA04 photodiode is ideal for measuring both pulsed and CW fiber light sources, by converting the optical power to an electrical current. The detector is housed in a TO-46 (mod) package with an anode, cathode, and case connection and mounted in an FC bulkhead connector. The photodiode anode produces a current, which is a function of the incident light power and the wavelength. The responsivity

 $\Re(\lambda)$, can be read from Figure 1 to estimate the amount of photocurrent to expect. This can be converted to a voltage by placing a load resistor (R_L) from the photodiode anode to the circuit ground. The output voltage is derived as:

$$V_o = P \times \Re \times R_I$$

The bandwidth, f_{BW} , and the rise time response, t_R , are determined from the diode capacitance, C_J , and the load resistance, R_L , as shown below. The diode capacitance can be lowered by placing a bias voltage from the photodiode cathode to the circuit ground. The effect of the bias voltage verse diode capacitance can be read from Figure 2.

$$f_{BW} = \frac{1}{(2\pi)R_L C_J}, t_R = \frac{0.35}{f_{BW}}$$

Specifications

Specification	Value	Specification	Value
Wavelength Range	800 - 1800 nm	Sensor Material	InGaAs Pin
Active Area	Ø100 μm (.0079 mm ²)	Max Ratings	
Rise/Fall Time (R_L =50 Ω , 5 V)	0.1 ns Min	Damage Threshold CW	70 mW
NEP, Typical (1550 nm)	$1.5 \text{ x } 10^{-15} \text{ W/Hz}^{\frac{1}{2}}$	Max Bias (Reverse) Voltage	25 V
Dark Current (5 V)	0.5 nA Typ. (1.0 nA Max)	Operating Temperature	-40 to 85 °C
Capacitance (5 V)	1.0 pF Typ. (1.2 pF Max)	Storage Temperature	-40 to 125 °C
Package	TO-46 (mod) w/ FC/PC	Reverse Current	10 mA
	Bulkhead Connector	Forward Current	10 mA

Recommended Circuit Diagram



* Case ground for PD with a third lead.





FGA04

Product Specification Sheet



Graphs



Capacitance vs Reverse Bias



Figures 1 and 2

Product Specification Sheet

THORLADS

Drawings



USA, Canada, & S. America Thorlabs, Inc. 435 Route 206 Newton, NJ 07860, USA Tel: 973-579-7227 Fax: 973-300-3600 www.thorlabs.com email: feedback@thorlabs.com

Europe Thorlabs GmbH Hans-Böckler-Str. 6 85221 Dachau, Germany Tel: +49-(0)8131-5956-0 Fax: +49-(0)8131-5956-9 www.thorlabs.com email: Europe@thorlabs.com UK and Ireland Thorlabs LTD. 1 Saint Thomas Place, Ely Cambridgeshire CB7 4EX, GB Tel: +44 (0)1353-654440 Fax: +44 (0)1353-654444 www.thorlabs.com email: sales.uk@thorlabs.com

Scandinavia Thorlabs Sweden AB Box 141 94 400 20 Göteborg, Sweden Tel: +46-31-733-30-00 Fax: +46-31-703-40-45 www.thorlabs.com email: scandinavia@thorlabs.com

Japan and Asia Thorlabs Japan Inc. 5-17-1, Ohtsuka Bunkyo-ku, Tokyo 112-0012, Japan Tel: +81-3-5979-8889 Fax: +81-3-5979-7285 www.thorlabs.jp email: sales@thorlabs.jp



3154-S01 Rev G – December 6, 2010 Specifications subject to change without notice.